

# Mobile Bioprocessing Containers 200L



## Engineering Excellence in Single-Use Bioprocessing...it's in our DNA

**CLARIPURE**<sup>®</sup> Bioprocessing containers for Clean room utilization. A range of plastic containers, manufactured in Germany, by AUER.



### Mobile Plastic Bioprocessing Containers:

External max. length (L)	810 mm
External max. width (W)	610 mm
External max. height (H)	969 mm
Internal min Length (Li)	738 mm
Internal min. Width (Wi)	537 mm
Internal max. height	658 mm
Individual Load Capacity	350 kg
Max. Static Load on bottom container	1500 kg
Max. Dynamic Load on bottom container	600 kg
Max. Stack Static	1+3
Max. Stack Dynamic	1+1
Collapsed unit height	530mm
Double stacked	1020 mm
Volume	200 L
Net weight kg	24,5 kg
Operational Temperature Range	-15 °C - 40 °C
BD position	variable