

Mobile Bioprocessing Containers 500L



Engineering Excellence in Single-Use Bioprocessing...it's in our DNA

CLARIPURE® Bioprocessing containers for Clean room utilization. A range of plastic containers, manufactured in Germany, by AUER.



Mobile Plastic Bioprocessing Containers:

External max. length (L)	1215 mm
External max. width (W)	805 mm
External max. height (H)	1086 mm
Internal min Length (Li)	1117 mm
Internal min. Width (Wi)	718 mm
Internal max. height	695 mm
Individual Load Capacity	600 kg
Max. Static Load on bottom container	4000 kg
Max. Dynamic Load on bottom container	1200 kg
Max. Stack Static	1+3
Max. Stack Dynamic	1+1
Collapsed unit height	575mm
Double stacked	1130 mm
Volume	500 L
Net weight kg	96 kg
Operational Temperature Range	-15 °C - 40 °C
BD position	Variable